



Application Note

AN 1013.00

Averaging in Color Measurement

“For most samples, a combination of optical and statistical averaging yields the best quantification of sample color.”

ABSTRACT

Averaging is a common colorimetric tool to “even out” the measurement highs and lows that are possible within a single sample or batch. The average measurement best represents the color of the sample or batch as a whole.

AVERAGING CAN TAKE TWO FORMS - OPTICAL AND STATISTICAL

Colorimetric instruments optically “average” all spatial detail within the field of view of the sensor to provide a color reading representing the average color within the view. It is recommended that the largest area of sample view possible be measured. The reason is that when you double (2x) the diameter of the sample view, the optical averaging is tripled (3x) based on the sample area. In addition to optical averaging for a single reading, multiple readings can be statistically averaged for a single measurement that best represents the color of the sample.

This application note further considers statistical averaging as a tool.

TYPES OF SAMPLES FOR AVERAGING

The ideal sample for yielding repeatable color measurements is completely opaque or transparent, flat, smooth and homogenous. However, as few samples meet this ideal, a variety of sample preparation and presentation techniques are employed to make color measurements as repeatable as possible. One technique for reporting the best result is to statistically average multiple readings from the same sample in different areas of the sample (measuring with replacement),

or average multiple readings of different samples from the same lot.

To illustrate, shown in Table 1 are five (5) individual readings made on a D25LT tristimulus colorimeter from a bag of pretzel pieces in a 5 inch diameter sample cup. The range and standard deviation provide an estimate of the measurement variation on this very non-uniform sample.

In Table 2, five (5) readings were averaged for each measurement, reducing the measurement variability. Between readings, the sample cup was emptied and refilled from the same bag. A single reading alone quantifies the approximate color. The average of the five readings best represents the overall color of the bag.

Examples of samples that benefit from averaging include non-homogenous translucent liquids with suspended particles, samples containing bubbles, scratches, and hazy areas, those with random color or texture distributions and directional samples such as corduroy fabric, ruled paper, and yarn. In each case, the samples should be read at least two to four times, with rotation and/or refills between each readings so that the sample is well presented to the instrument.

The statistical average of the multiple readings, using the largest area of sample view possible, is reported as best quantifying the color of the sample.

TABLE 1. 5 INDIVIDUAL READINGS OF PRETZEL BITES

Single Reading per Measurement	CIE C/2		
	L*	a*	b*
1	42.79	8.57	38.57
2	42.25	6.88	38.43
3	42.88	8.03	33.83
4	41.03	8.98	33.09
5	43.11	9.62	35.27
Range	2.08	2.74	5.48
Standard Deviation	0.83	1.04	2.55

TABLE 2. STATISTICAL AVERAGE OF 5 READINGS OF PRETZEL BITES WITH REPLACEMENT

Average of 5 Readings per Measurement	CIE C/2		
	L*	a*	b*
1	41.55	8.75	37.11
2	41.91	9.33	34.70
3	42.12	9.33	33.27
4	41.96	9.08	33.31
5	42.28	8.76	33.28
Range	0.73	0.58	3.84
Standard Deviation	0.27	0.29	1.67



Figure 1. Glass sample cup with pretzel bites ready for color measurement on a HunterLab D25LT with a large 89 mm diameter area of sample view.



Figure 2. This is the instrument view through the bottom of the dish. The instrument optically averages all colors within the field of sample view.

DEFINE THE MEASUREMENT METHOD

Six parameters are important for any good color measurement method and optical and statistical averaging is a key part of that method. Changing any of the following parameters will affect the color values.

1. **Color Scale** - Select the color or color difference scale. If unsure, use CIE L^* , a^* , b^* or Hunter L , a , b .
2. **CIE Illuminant** - Select the illuminant. If unsure, use CIE D65. The D25LT is fixed for the CIE illuminant C representing average daylight.
3. **CIE Standard Observer** - Select the observer. If unsure, use the 1964 10 degree Standard Observer. The D25LT is fixed to the 1931 2 degree Standard Observer.
4. **Instrument Geometry** - This parameter is fixed by the choice of instrument model which can be a directional $45^\circ/0^\circ$ or diffuse $d/8^\circ$ sphere. Typically an instrument with a directional $45^\circ/0^\circ$ like the D25LT is used to measure the very non-uniform samples.
5. **Sample Preparation** - Sample preparation must be consistent with techniques that support making the sample more uniform leading to a more repeatable reading.
6. **Sample Presentation** - Presentation techniques use optical (largest area of sample view) and statistical (multiple readings with replacement per measurement) to ensure a repeatable and representative color measurement of the sample.





Figure 3. Grooved Kraft board showing directionality effect with 90° rotation.

CONCLUSION

The variation between measurements can be lowered by averaging multiple readings per measurement with replacement, using the largest area of view possible with each reading. It is important to document the measurement methodology for the purpose of effectively replicating and communicating color values.



*More Information about
Color Measurement on our
HunterLab Blog*

measuretruecolor.com

ABOUT HUNTERLAB

HunterLab, the first name in color measurement, provides ruggedly dependable, consistently accurate, and cost effective color measurement solutions. With over 6 decades of experience in more than 65 countries, HunterLab applies leading edge technology to measure and communicate color simply and effectively. The company offers both diffuse/8° and a complete line of true 45°/0° optical geometry instruments in portable, bench-top and production in-line configurations. HunterLab, the world's true measure of color.

© Hunterlab 2012

Hunter Associates Laboratory Inc.,
11491 Sunset Hills Road, Reston, VA 20190-5280 USA
helpdesk@hunterlab.com
www.hunterlab.com

