

Fundamentals of Color and Appearance

Module 9

Haze, Turbidity, Opalescence

Color Science Educational Series



www.hunterlab.com | support.hunterlab.com



In manufacturing environments where visual appearance matters, terms like clarity and cloudiness are commonly used to describe how well light passes through a material or liquid.



Haze, Turbidity, Opalescence



While “clarity” is not a standardized or directly measurable scientific parameter, it serves as a practical and intuitive way to describe the visual impression of transparency.

Conversely, “cloudiness” reflects the presence of visible or sub-visible light scatter, often perceived as haze or milkiness.



Fundamentals of Color and Appearance

Module 9:

Haze, Turbidity, Opalescence



To move from subjective descriptions to objective, quantifiable measurement, we rely on three distinct but related scientific optical properties.




Haze, Turbidity, Opalescence



These are **Haze**, **Turbidity**, and **Opalescence**. At the heart of all three properties is the principle of light scatter.

Let's begin with a definition and everyday example for each property.





Let's start with **Turbidity**.

Imagine standing by a clear stream on a calm day. The water is transparent. You can see the pebbles on the bottom, the fish swimming, even the shadows they cast.





Now picture that same stream after a heavy rain.

The water is muddy and murky from Suspended soil, silt, and debris that scatter and absorb light. This is turbidity in action.





This is an important measurement parameter for many liquids including food and beverage manufacturing.





Now let's talk about **Haze**.

Imagine gazing out at a mountain range on a clear day. You can see the contours of each mountain ridge and even the shadows on the slopes. Every detail is sharp and vivid.



Haze, Turbidity, Opalescence



Now imagine the same view on a hot, humid summer day or during wildfire season.

A thin veil of atmospheric haze made of dust, moisture, or smoke fills the air. The mountains are still visible, but their edges blur, colors fade, and everything looks washed out.





This is how haze affects transparent materials. Scattered light blurs the view, reducing contrast and clarity. This is an important measurement parameter for applications like food and beverage packaging.





And Haze is an important measurement parameter for ophthalmic lenses used in eyeglasses.





Finally, let's look at Opalescence.

Imagine walking by a lake that appears sharp, distinct and bright blue.



Haze, Turbidity, Opalescence



Now imagine walking by a lake that looks like it has a soft glow, maybe a gentle teal or white shimmer. This is opalescence.

It's not the full-on cloudiness of turbidity, It's a slight, pearly scatter of light, often caused by microscopic particles. It's most noticeable when light passes through at the right angle.





In pharmaceutical manufacturing, especially for injectable drugs and biologics, this delicate visual disturbance can signal a critical issue.



Haze, Turbidity, Opalescence



Though all three properties relate to light scattering, they are quantified using different scientific measurement methods.





Turbidity is measured using nephelometers or turbidimeters.

The results, in NTU or nephelometric turbidity units, help determine whether a sample is acceptably clear, especially critical for water, clear beverages, and some injectable solutions.





Haze is used to evaluate the amount of light scatter in transparent solids using ASTM method D1003, which outlines two optional procedures:



Haze, Turbidity, Opalescence



ASTM D1003:

- **Procedure A** uses a dedicated haze meter and an integrating sphere to compare total transmitted light with scattered light beyond a 2.5° angle.
- **Procedure B** allows use of a spectrophotometer with an integrating sphere—ideal for labs already measuring color and transmittance.

In both cases, the result is expressed as a percentage of scattered light, with lower values indicating clearer, higher-quality materials.






It is important when measuring haze to understand the difference between **Total Transmission**, or **TTRAN**, and **Regular Transmission**, or **RTRAN** for accurate measurements.

Total Transmission is the total amount of light that passes through a transparent sample. It includes both direct and diffusely scattered light at all angles.

Regular Transmission measures only the direct, non-scattered light that passes through the sample in a straight line.


- **TTRAN (Total Transmission)** is the total amount of light that passes through a transparent sample.
 - It includes both direct (unscattered) light and diffusely scattered light at all angles.
- **RTRAN (Regular Transmission)** measures only the direct, unscattered light that passes through the sample in a straight line (within a very narrow angle—typically $\pm 2.5^\circ$ of the incident beam).




Haze, Turbidity, Opalescence



Total Transmission (TTRAN) is measured by placing the sample flush against the sphere port. This captures the total amount of light that passes through a transparent sample and includes both direct light and diffusely scattered light at all angles.



Regular Transmission (RTRAN) is measured by placing the sample anywhere away from the sphere port, capturing only the direct light that passes through the sample in a straight line.





Separating **direct** vs. **scattered** light transmission allows precise quantification of haze, critical for evaluating optical clarity in plastics, films, lenses, and other transparent materials.

A screenshot of a software interface displaying a table of color and haze data. The table has columns for Name, Date, Time, L*, a*, b*, and Haze % [C/2]. The data rows include a Standard, Tolerances, and three samples (Sample 3, Sample 2, and Sample 1). A green "Measure" button is visible at the bottom right of the interface.

Name	Date	Time	L*	a*	b*	Haze % [C/2]
Standard (Ad Hoc / Working)			95.48	0.09	0.44	32.68
Tolerances	--	--	--	--	--	--
▶ Sample 3	09/02/2025	15:55	95.49	0.09	0.43	32.66
▶ Sample 2	09/02/2025	15:55	95.48	0.09	0.44	32.67
▶ Sample 1	09/02/2025	15:54	95.48	0.09	0.44	32.68



Lastly, **Opalescence**, typically assessed in pharmaceutical quality control using visual comparison, is guided by pharmacopeial standards like the European Pharmacopoeia which defines reference suspensions at known NTU levels.





For greater accuracy and precision however, instrumental nephelometry is used, as defined in USP method 855 or European Pharmacopoeia method 2.2.1.

Instruments are calibrated using liquid Formazin suspensions, just like turbidity measurements, but with greater sensitivity for very low turbidity levels.

Opalescence is typically assessed in pharmaceutical quality control using Visual comparison, guided by pharmacopeial standards



Haze, Turbidity, Opalescence



Understanding, measuring and controlling haze, turbidity and opalescence isn't just about aesthetics, it's about safety, performance and brand trust.





In consumer packaging, haze affects brand perception and purchasing decisions.





In beverages, turbidity can signal quality or spoilage.





And in injectables, opalescence can mean the difference between a stable, safe drug and one that must be rejected.





Each property acts as a quality control checkpoint, whether in material production, formulation, or final product inspection.





Thank you for joining us on this journey through color science.

Be sure to watch Module 10, **Color verses Appearance** next, where we will explore how a materials surface characteristics like texture and gloss not only influence how we perceive color but also present real challenges in accurately measuring color..

And be sure to visit hunterlab.com to learn more about how our solutions can help you achieve color confidence, every time, or to schedule a consultation with one of our color experts.

Fundamentals of Color and Appearance

Module 10

Color vs. Appearance

Color Science Educational Series



www.hunterlab.com | support.hunterlab.com

**Schedule a free color consultation
with our experts**



www.hunterlab.com