



# Enhancing Cracker Manufacturing with Spectrophotometric Color Measurement

## Introduction

Crackers are a staple in the global snack food market, valued for their crisp texture, light golden appearance, and uniformity. Consumers often buy crackers in bulk or multipacks, where visual consistency across every cracker is essential for maintaining confidence in the brand. Even subtle variations in surface color, whether from uneven baking, excess salt application, or ingredient inconsistencies—can signal poor quality to consumers.

Spectrophotometric color measurement offers an objective, repeatable way to monitor cracker appearance throughout production. By quantifying color, manufacturers can reduce waste, standardize quality, and strengthen brand trust. This white paper explores why color is a critical quality parameter in crackers, the challenges of visual inspection, and how **HunterLab's Aeros spectrophotometer** provides the most accurate, efficient solution.

## Importance of Color in Crackers

- **Indicator of Bake Quality** - Golden hues communicate crispness and flavor; pale color suggests underbaking, while dark brown indicates overbaking.
- **Ingredient and Topping Consistency** - Crackers are often topped with salt, seeds, or seasonings. Excess or uneven application can change surface reflectance and consumer perception.
- **Freshness Cue** - Uniform lightness reassures consumers the crackers are properly baked and not stale.



- **Brand Identity** - Each cracker style (e.g., salted, whole wheat, multigrain) has a defined appearance that consumers expect.

## Challenges in Cracker Appearance Control

- **Heterogeneous Surfaces** - Salt crystals or toppings scatter light unpredictably, making visual judgments unreliable.
- **Oven Variability** - Slight differences in dwell time, temperature, or airflow across the oven belt affect browning.
- **Human Subjectivity** - Visual QC depends on lighting and perception, which vary between inspectors.

## Why Instrumental Color Measurement is Essential

Spectrophotometric measurement eliminates subjectivity by reporting precise *CIELAB values* ( $L^*$ ,  $a^*$ ,  $b^*$ ) under controlled conditions.

- $L^*$  (lightness) is especially useful for crackers, as it correlates with both surface toppings and bake level.
- Multiple readings averaged across a tray or batch ensure representative results.
- Results can be tied to clear pass/fail tolerances, supporting global quality consistency.

## Recommended Solution - HunterLab Aeros

The HunterLab Aeros is ideal for crackers:

- **Non-Contact, Large Area Measurement** - Measures trays or piles of crackers without grinding, compressing, or cleaning.



- **27.5 in<sup>2</sup> Capture with 35 Readings in 5 Seconds** - Averages across multiple crackers for highly representative results.
- **Auto Height Sensing** - Adapts to varying sample depths for consistent readings.
- **Integrated Touchscreen and Software** - Provides immediate pass/fail analysis and data storage.

#### Benefits for Cracker Manufacturers:

- Detects baking deviations early.
- Identifies topping inconsistencies (excess salt, seeds, etc.).
- Reduces consumer complaints and costly rework.
- Strengthens brand confidence through objective QC.

## Hypothetical Case Studies

### Case Study 1 - Excessive Salt Detected by L\* Value

#### Background

A premium cracker brand applies coarse salt as a topping. Consumers expect a light golden cracker with evenly distributed white salt crystals.

#### Challenge

During production, the **Aeros** spectrophotometer detected significantly higher **L\*** values compared to the standard. While the dough and bake were consistent, the crackers appeared lighter overall.

#### Investigation

Analysis revealed that the salt dispenser was miscalibrated, applying too much salt per cracker. The reflective salt crystals increased surface brightness, driving up the average **L\*** reading.



## Corrective Action

- The salt feeder was recalibrated to restore correct topping levels.
- Aeros checkpoints were placed post-topping to continuously monitor L\* values before baking.

## Outcome

The manufacturer prevented excessive-salt crackers from reaching packaging, reducing potential customer complaints of “over-salted” product and ensuring consistent appearance.

## Key Takeaway

Aeros can detect not only bake differences but also topping-related inconsistencies, making L\* a powerful early-warning metric for multiple quality attributes.

## Case Study 2 – Over Doneness Detected by L\* Value

### Background

Another batch of the same crackers showed an unexpected drop in L\* values, indicating darker product compared to the golden-brown standard.

### Challenge

Visually, the crackers appeared overbaked, with some edges nearing a burnt tone. If packaged, these could be perceived as stale or poor quality.

### Investigation

Aeros readings confirmed the crackers were consistently outside the  $\Delta E$  tolerance. Process review revealed a drift in oven temperature settings, leading to excessive browning.

## Corrective Action



- Oven controls were recalibrated to maintain proper bake conditions.
- An Aeros checkpoint was added immediately to the post-oven to catch over-bake events early.

## Outcome

Only a limited quantity was affected, saving the manufacturer from widespread rework or consumer rejection. Once corrected, subsequent batches aligned perfectly with the target golden-brown  $L^*$  profile.

## Key Takeaway

For crackers,  $L^*$  value provides a sensitive, objective measure of proper bake level. Aeros enables real-time detection of over-doneness, preventing brand damage and product waste.

## Conclusion

Cracker manufacturing depends on consistent appearance to communicate freshness, quality, and brand identity. Both surface toppings and bake level influence consumer perception, and visual QC alone cannot reliably capture these variations.

By implementing **HunterLab's Aeros spectrophotometer**, manufacturers gain:

- Objective monitoring of golden-brown targets ( $L^*$ ,  $a^*$ ,  $b^*$ ).
- Fast, representative results across heterogeneous cracker surfaces.
- Early detection of both topping inconsistencies and baking deviations.

With Aeros, cracker producers can deliver every box with confidence—ensuring that crackers look as good as they taste, reinforcing consumer trust, and protecting brand reputation.