

## HunterLab Sample Trays for D25 NC and D25 NC FM



*D02-1016-490 Assembly, Sample Tray, 12" Diameter x 1" Depth  
Shipped standard with the D25 NC*



*D02-1017-314 Sample Tray, Nested 12" x 1" & 6" x 1"*

### **Anodized Sample Trays**

HunterLab offers a selection of dark anodized sample trays. These trays come in a variety of diameters and sidewall heights to best fit your sample measuring needs.

Optional nested tray assemblies are available for purchase. This assembly consists of a smaller tray inside a larger diameter tray. This option reduces the amount of sample required in the space between the two trays. It is not possible to nest individually sold pans.

Note: HunterLab has previously offered brushed aluminum sample trays. We believe these new trays will provide greater performance through longevity and ease of maintenance. The different color of trays will not affect your sample readings.

### **Anodized Standard Sample Trays**

The D25 NC will be shipped with two standard trays (No exchanges available):

D02-1016-971	Assy, Sample Tray, 6" Dia x 1" Depth
D02-1016-490	Assy, Sample Tray, 12" Dia x 1" Depth

As always, no trays are shipped as standard with the D25 NC FM instrument. You may order the L02-1014-741 Sample Handling Package for Chips as an accessory for the D25 NC FM.

The following anodized trays are available to be purchased separately through HunterLab:

D02-1016-971	Assy, Sample Tray, 6" Dia x 1" Depth
D02-1017-308	Assy, Sample Tray, 8" Dia x 1" Depth
D02-1017-309	Assy, Sample Tray, 10" Dia x 1" Depth
D02-1016-490	Assy, Sample Tray, 12" Dia x 1" Depth
D02-1017-310	Assy, Sample Tray, 6" Dia x 2" Depth
D02-1017-311	Assy, Sample Tray, 8" Dia x 2" Depth
D02-1017-312	Assy, Sample Tray, 10" Dia x 2" Depth
D02-1017-313	Assy, Sample Tray, 12" Dia x 2" Depth
D02-1017-314	Sample Tray, Nested 12" x 1" & 6" x 1"
D02-1017-315	Sample Tray, Nested 12" x 2" & 6" x 2"

**All part information and pricing may be found on the D25 NC Price List.**