

ACCESSORIES

ColorFlex[®] L2

HunterLab accessories are designed to enhance the functionality and precision of your color measurement instruments. From sample handling tools to protective components, these accessories ensure accurate results and seamless operation.

What's in the Box

All you need to get accurate color measurement in minutes!



- ColorFlex L2 Instrument • Standards Case • Diagnostic Check Tile (Green Tile) • Instrument Standard (White Tile)
- Reflectance Black Glass • Port Plate 25 mm (1 in) • Lens Wipes Power Supply • Quick Start Guide
- Stylus with Holder • Cleaning Cloth • Flash Drive • Remote Access Support Flyer

OPTIONAL ACCESSORIES



GLASS SAMPLE CUP 64mm (2.5 inch)

Optically clear glass cup for sample presentation of liquids, powders, granules and pellets.

04-7209-00



RING AND DISK SET

For use in the glass sample cup for sample presentation of liquids, powders, granules and pellets.

02-4579-00

ACCESSORIES



OPTIONAL ACCESSORIES



SAMPLE CUP OPAQUE COVER

Provides a light trap to exclude the interference of external light on the sample cup.

04-4000-00



OPEN PORT INSERT for Sample Cup

The open port insert centers the 64 mm (2.5") sample cup at the measuring port, with a drop-down lip for lens focus on the glass window's bottom.

04-6622-00



OPEN PORT INSERT 13 mm (0.5")

The open port insert is a field-replaceable, snap-in component that reduces the instrument's 25 mm viewing area to a smaller 13 mm diameter.

04-6623-05



OPEN PORT INSERT 19 mm (0.75")

The open port insert is a field-replaceable, snap-in component that reduces the instrument's 25 mm viewing area to 19 mm.

04-6623-08



OPEN PORT INSERT 25 mm (1")

The open port insert provides a field-replaceable, snap-in 25 mm (1") port insert that is open and masks the sample to the 25 mm viewing area of the instrument.

04-6623-09



OPEN PORT INSERT 25 mm (1") with Glass

The port insert is a field-replaceable, snap-in 25 mm (1") component with clear borosilicate glass, masking the sample to a 25 mm viewing area while protecting the instrument from particles.

02-6624-01



OPEN PORT INSERT 32 mm (1.25") with Glass

The port insert is a field-replaceable, snap-in 32 mm (1.25") component with borosilicate glass, masking the sample to a 32 mm viewing area and protecting the instrument from particles.

02-6624-05

For more information, please contact HunterLab at 703-471-6870, sales@hunterlab.com or visit www.hunterlab.com

HunterLab and ColorFlex are Trademarks of Hunter Associates Laboratory, Inc.

The accessories listed are compatible with ColorFlex L2 as specified. Availability, specifications, and features are subject to change without notice.