



Compact
Powerful
Precise



The world's true measure of color

HunterLab



ColorFlex L2

*Level UP your Color Confidence with the new
ColorFlex L2 Spectrophotometer.*




Explore the Next-Gen ColorFlex®...L2

As simple as you want, as powerful as you need.

Elevate your color accuracy and user experience like never before. Designed to streamline and expedite your quality control processes, ColorFlex L2 Spectrophotometer gives you unsurpassed confidence and trust in your outcomes. Compact and ergonomic build requires minimal space while offering enhanced durability and reliability.



Explore all the additional features and benefits below.



45°a:0° Annular Illuminated Optical Geometry

Superior Color Accuracy and Consistency: Enhancing the efficiency of Quality Control workflows and delivery of superior quality products repeatedly.



Compact with Sealed Spillproof Case

Space-Saving and Easy Maintenance: Occupies minimal space and is designed to withstand the demands of busy lab environments.

Sliced Cheese

SPICES

Puree

TOMATO KETCHUP

Yogurt

Tomato Paste

flour

mustard

barium sulfate

Baking Powders

Citrus Juice

Herbs

CEMENT




Powerful Standalone Functionality

No computer required: With an onboard computer and huge storage capacity, ColorFlex L2 can function as an independent color workstation. Connect monitor, keyboard and mouse if desired.



Modern Tablet User Interface with Touchscreen & Color Display

Durable and Efficient: Facilitates accurate color measurement in seconds, with a robust industrial glass design suitable for harsh environments.




New & Improved Color Software

Essentials 2.0: Take precise and reliable color measurement in seconds with next-gen user-friendliness and seamlessly analyze multiple tristimulus and spectral data sets. It's that powerful.

Self-Guided and Intuitive Out-of-Box Startup

Time-Saving and User-Friendly: Ready to use straight out of the box with an integrated wizard for quick training, reducing operator time and easing concerns related to staff training and turnover.



Powdered Sugar
tomato juice
RAW GRAINS

BAKING CHIPS
Magnesium Sulfate
Starch
powdered cheese


Powdered Functional Beverages
shortening
Orange Juice
Tomato Soup

WHITE SUGAR



Built-in Camera for Sample Viewing and Screen Recording

Precision and Communication: Ensures accurate sample positioning to eliminate measurement discrepancies and allows easy sharing of recordings for swift quality control decisions.



Built-in Advanced Diagnostics

On-Board Diagnostic Health Check: Ensures the HunterLab ColorFlex L2 is always operating at peak performance, continuing to build on the legendary history of reliability of the ColorFlex family of instruments.

Advanced Data Export and Communication

Full suite of data communication options: Features USB ports, HDMI, and support for mouse and keyboard, enabling easy data export/printing.



Connectivity with LIMS and SPC systems

Simplifies QC workflows: Aids in trend analysis and workflow simplification across manufacturing sites.



- Tea Drinks (powder)
- Calcium Carbonate
- POWDERED SPORTS DRINKS
- Fertilizers
- BABY FOOD
- Tomato Sauce
- Brown Sugar
- Coffee - Green beans, Roasted & Ground
- SEASONINGS
- Concrete
- Cocoa Powder



ColorFlex L2 Spectrophotometer helps you reduce wastage of your end products, contribute towards your sustainability goals, and boost brand image and bottom-line revenue.

Welcome **the next-gen ColorFlex**, without any interruptions!



From **OLD...**



...to **NEW**

We have upgraded the most relied upon $45^{\circ}/0^{\circ}$ spectrophotometer in the world, ColorFlex EZ, to the new ColorFlex L2. Offering you an even better user experience without any interruptions! You can easily transfer your ColorFlex EZ standards and data to the ColorFlex L2. Same Xenon light source and $45^{\circ}a:0^{\circ}$ Annular Illuminated Optics as the ColorFlex EZ provide seamless transition from old to new.

Despite its small size, requiring very little benchtop space, the ColorFlex L2 is powerful and designed to meet your specific color measurement needs.

The HunterLab Advantage

The ColorFlex L2 is backed by over 70 years of quality innovation and experience from HunterLab, the world's most trusted color quality experts.

With an unmatched reputation for delivering the right solution, HunterLab tailors products and technologies for every color measurement need and budget offering the broadest range of color measurement solutions in the industry.

Our industry-leading clients trust HunterLab spectrophotometers to benefit from the functionality and accuracy needed for color quality. For more information or to schedule a demo of ColorFlex L2, visit <https://www.hunterlab.com/en/contact/>





CONTACT US

We've worked with thousands of companies, but you'll feel like our one and only. If color accuracy is crucial to your business, you've come to the right place—with the right people. Let's get in touch. Ask a question, Request a quote or just say hello!

OUR OFFICES

HunterLab

11491 Sunset Hills Road
Reston, VA 20190-5280
Website: www.hunterlab.com
Email: sales@hunterlab.com
Ph: 703-471-6870
Fax: 703-471-4237

Corporate Headquarters

HunterLab Europe

Christian Jansen, Managing Director
Dr.-August-Einsele-Ring 15
82418 Murnau, Germany
Website: www.hunterlab.de
Email: Christian.Jansen@hunterlab.com
Ph: +49 8841 9464

European Office

HunterLab Asia-Pacific

Leo Cheung, Regional Manager
Flat E, 13F Leahander Centre
28 Wang Wo Tsai Street Tsuen Wan
New Territories Hong Kong
Website: www.hunterlabcolor.cn
Email: sales-asiapacific@hunterlab.com
Ph: +852-2406 1723

Asia-Pacific Office



For more information, please contact HunterLab at 703-471-6870, sales@hunterlab.com or visit www.hunterlab.com