

EasyMatch QC-ER version 4.86 and higher Master List of Event Log codes

Event ID	Description	Category
1001	Modified Application Preferences	Application Options
1002	Modified Naming Conventions	Application Options
1003	Modified Default Tolerances	Application Options
1004	Accessed Group Privileges	Application Options
1005	Modified Server settings	Application Options
Authentication		
3001	Logged in	Authentication
3002	Login Failed	Authentication
3003	Login tried and cancelled	Authentication
3004	Logged Off	Authentication
Operations		
4003	Sample Read performed	General
4004	Description: Standard Read performed	General
4005	Saved Job file <path of the Job file>	General
4006	Closed Job file <path of the Job file>	General
4007	Job Saved to database : <Job Name>	General
4009	a) Series recalled from database:<Details of Series> b) Standard Recalled from database:<Details of Standard> c)Sample recalled from database:<Details of Sample>	General
4010	a) Series deleted from database b) Standard or Sample deleted from database	General