



Keyboard (A13-1014-294)

Instruments: ColorFlex, MiniScan

Purpose: The optional flexible keyboard connects to the instrument's USB port (via Standard-A to Mini-A adapter for the MiniScan (A21-1013-859)) to allow easy entry of any alphanumeric information (such as setup names) required.



The **UP, DOWN, LEFT, AND RIGHT** arrow keys on the keyboard correspond to the same buttons on the button pad. The **ENTER** key on the keyboard emulates the center (Go) button on the button pad. Numbers and letters may be typed on the keyboard as usual.

Installation: The keyboard is automatically detected by the instrument once connected.