

Ideal Sample For Color Measurement

 Opaque or transparent

 Solid

 Flat

 Smooth

 Uniform

Sample Preparation and Presentation

- ☑ Choose samples that are representative of the product.
- ☑ Prepare the sample in a way to best approximate the ideal sample characteristics.
- ☑ Prepare samples in the same way each time.
- ☑ Present the samples to the instrument in a repeatable manner.
- ☑ Take multiple readings of the sample for an average measurement.



Thank You



If you are connected to the internet click:



For more information about sample preparation and presentation



For product information



visit our website at www.hunterlab.com

or

contact us at info@hunterlab.com

