

# Addendum to MiniScan EZ User's Manual

---

## Description

CMR-3114 consists of the following extra parts for the MiniScan EZ :

1. Ball transfer roller feet (2) to facilitate movement of the MSEZ across a surface.
2. A bottom bracket to attach the roller feet.
3. A strain relief clamp which attaches to a nut plate on the housing by means of a screw.
4. Rubber dust cover for the connector opening.



**Metal bracket with ball transfer roller feet.**



**Strain relief clamp with rubber dust cover.**

## Installation

All elements of this CMR are factory installed.

## Operation

Operate the MiniScan normally.